

PLACEBO SUCKABLE TABLETS WITH NEOSORB® P60W SORBITOL AND XYLISORB® P300 XYLITOL

FORMULA

	% (by weight)	mg / tablet
NEOSORB® P60W sorbitol	76.73	1 381.1
XYLISORB® P300 xylitol	20.00	360.0
Magnesium stearate	1.20	21.6
Cherry flavour, Haarmann Reiner	1.00	18.0
Anhydrous citric acid	1.00	18.0
Red Ponceau lake, Sensient Colors	0.07	1.3
		1 800.0

METHOD

- In a mixer, mix NEOSORB® P60W sorbitol, XYLISORB® P300 xylitol, flavour, colouring agent and anhydrous citric acid for 20 minutes.
- Add magnesium stearate and mix for 5 minutes.
- Compress the powder in a single punch press fitted with 18mm flat bevel edge punches.

COMMENT

- Weight 1800 mg
- Thickness 5.02 mm
- Density 1.412
- Erweka hardness 464.1 N